DCHHS Submission Instructions for Novel Coronavirus (COVID-19) PCR Testing

Interim Guidance updated as of March 13, 2020

Criteria for Novel Coronavirus (COVID-19) Testing Approval

 Per the Texas Department of State Health Services (DSHS), as of 3/4/20, all public health laboratories in Texas will continue to use the CDC PUI criteria¹ below posted on February 27, 2020 to prioritize testing of persons at higher risk for COVID-19: <u>https://emergency.cdc.gov/han/2020/han00428.asp</u>

| Texas DSHS Criteria for Testing (as of 3/11/20) | | |
|---|-----|---|
| Clinical Features | & | Epidemiologic Risk |
| Fever or signs/symptoms of lower respiratory illness (e.g., cough or shortness of breath) | AND | Any person, including health care workers, who has had close contact with a laboratory- confirmed 2019-nCoV patient within 14 days of symptom onset |
| Fever and signs/symptoms of a lower respiratory illness (e.g., cough or shortness of breath) | AND | A history of travel from <u>affected geographic</u> <u>areas</u> within 14 days of symptom onset OR An individual(s) with risk factors that put them at higher risk or poor outcomes |
| Fever and signs/symptoms of a lower respiratory illness (e.g., cough or shortness of breath) requiring hospitalization | AND | No source of exposure has been identified and other causes of respiratory illness have been ruled out (e.g., influenza) |

- To discuss testing of patients who reside in Dallas County, clinicians should contact the DCHHS Epidemiology division at 214-819-2004 during business hours or (877) 605-2660 after hours and complete the required <u>DCHHS Novel Coronavirus Test Request Form</u>. For patients who reside outside of Dallas County, submitter <u>must obtain approval</u> for testing from the <u>appropriate county/state health</u> department³ prior to submission to DCHHS.
- The DCHHS LRN Laboratory can only accept specimens for novel coronavirus PCR testing from residents of the following counties comprising its service area: **Collin, Dallas, Ellis, Fannin, Grayson, Henderson, Hunt, Kaufman, Navarro, Rains, Rockwall, and Van Zandt**.

Specimen Collection Instructions:

For initial diagnostic testing for COVID-19, the CDC recommends collecting the following specimen⁴:

Upper Respiratory Tract:

- Nasopharyngeal swab (NP swab)
 - Use only a synthetic fiber swab with aplastic shaft. Do not use a calcium alginate swab or swab with a wooden shaft, as they may contain substances that inactivate some viruses and inhibit PCR testing. Place swab immediately into a sterile tube containing 2-3 ml of viral transport media. Refrigerate specimen at 2-8°C and ship on ice packs.
 - <u>Nasopharyngeal swab</u>: Insert a swab into the nostril parallel to the palate. Leave the swab in place for a few seconds to absorb secretions.

[Whole blood, serum, urine, and stool will NOT be accepted.]



Criteria for Specimen Handling

- Specimen <u>MUST</u> be labeled with the following information:
 - 1. Patient's Full Name
 - 2. Patient's Medical Record Number (MRN)
 - 3. Patient's date of birth (DOB)
 - 4. Specimen Type (e.g. NP)
 - 5. Date of Specimen Collection
- A specimen that is not appropriately labeled will <u>NOT</u> be accepted for testing.
- Specimen should be double-bagged in a biohazard bag (one bag per patient) and must be
 accompanied by a completed printed copy of the <u>DCHHS Novel Coronavirus Test Request</u>
 <u>Form</u> placed in the outer pouch of the biohazard bag. Multiple patient samples (individual
 bags for each patient) may be batched together in a single bag for delivery.

Specimen Transport

- For patients who reside in Dallas County, submitter <u>MUST obtain approval from the DCHHS</u> <u>Epidemiology division</u> before a specimen can be sent to DCHHS for testing. Please consult with DCHHS Epidemiology at 214-819-2004 during business hours or (877) 605-2660 after hours. In addition to a paper copy of the completed <u>DCHHS Novel Coronavirus Test Request Form</u> in the specimen shipment, please also fax one copy to 214-819-1933.
- For patients who **reside outside of Dallas County**, submitter <u>must obtain approval</u> for PUI reporting and testing from the <u>appropriate respective county/state health department</u>³ prior to submission to the DCHHS LRN Laboratory.
- All healthcare facilities must arrange for transport of specimen from their facility to DCHHS; DCHHS will not pick up specimens from any submitters.
- Refrigerate specimen at 2-8°C and ship on cold packs within 24 hours of collection. After 24 hours of collection, specimen should be frozen at -70°C and shipped on dry ice.
- Follow packing and shipping instructions for Category B, Biological Substances (see https://www.iata.org/contentassets/b08040a138dc4442a4f066e6fb99fe2a/dgr-61-en-pi650.pdf for instructions).
- Ship according to IATA (International Air Transport Association) guidelines using overnight courier to arrive during business hours Monday–Friday, 8:00 am 4:30 pm.
- Specimens delivered by same-day courier services between Monday–Friday, 8:00 am 4:30 pm will also be accepted.

DCHHS LRN Shipping and Delivery Address:

Dallas County Health and Human Services Laboratory 2377 N. Stemmons Fwy., Basement Ste. 003 Dallas, TX 75207 Laboratory: 972-692-2762 Epidemiology: 214-819-2004 After-hours: 877-605-2660